

USA BlueBook BRAND PRODUCTS

YOUR BEST VALUE!



USA BlueBook Electrode Solutions

A Electrode Storage Solutions

To extend your pH electrode's life and ensure it's ready to go when needed, store it in our electrode storage solution between uses.

DESCRIPTION	STOCK #	EACH
Electrode Storage Solution, 500 mL	40450*	\$
Electrode Storage Solution, 1 L	201441*	

* Not for use with Hach IntelliCAL electrodes or PHC28101 (stock # 37772).

B Electrode Cleaning Solutions

Convenient premixed cleaning solution is ideal for removing crystallized deposits on your pH electrode bulb or reference junction. For clogged reference junctions, soak electrode in dilute, silver-free KCl solution at 70°C for 20 minutes. (Mix 3 mL of stock # 40394 with 100 mL distilled water for dilute solution.)

DESCRIPTION	STOCK #	EACH
Electrode Cleaning Solution, 500 mL	40453	\$
4M KCl Solution, 60 mL	40394	

C Electrode Fill Solutions

These solutions are compatible with many popular brand electrodes. Listed below are some equivalent brand electrode solutions. We recommend you replace, not just refill, your solution every 2 to 4 weeks to extend the life of your electrode. Sold in 500 mL poly-containers.

DESCRIPTION	EQUIVALENT BRANDS	STOCK #	EACH
Potassium Chloride, 3 M	Coming 477006 Hach SensiON+ Thermo Orion 810007 WTW SenTix	29731	\$
Potassium Chloride, 4 M (Saturated)	USABlueBook (Double-Junction) PlantPRO 72368 Oakton 35803-74	29732	
Potassium Chloride, 4 M, Sat. with Silver Chloride	USABlueBook (Single-Junction) Oakton 35803-73 Thermo Orion 900011	29733	

Electrode Fill Solutions

- For longest life, refill your electrode every 2 to 4 weeks



For Hach sensiON Electrodes

Use with: 5191000, 5191011, 5191500, 5191022

MFR/REFILL #	SIZE	STOCK #	EACH
2546902	2 Cartridges	44554	\$

For Hach sensiON+ Electrodes

MFR/REFILL #	SIZE	STOCK #	EACH
LZW9509.99	50 mL	204771	\$
LZW9510.99	125 mL	204772	

For Hach IntelliCAL Electrodes (3.0M KCl sat. with AgCl)

Use with: PHC30101, PHC30103

MFR/REFILL #	SIZE	STOCK #	EACH
2841700	28 mL	39171	\$

For Hach IntelliCAL Electrodes (2.44M KCl)

Use with: PHC28101

MFR/REFILL #	SIZE	STOCK #	EACH
2965029	59 mL	39175	\$

For Orion Electrodes (3M KCl)

Use with: 80-03, 80-05, 80-05U, 81-02, 81-02U, 81-03, 81-03U, 81-04, 81-04U, 81-15, 81-15U, 81-35, 81-35U, 81-55, 81-56, 81-56U, 81-57, 81-62, 81-63, 81-65, 81-66, 81-72, 81-75, 82-02, 82-03, 82-35, 82-56, 82-63, 83-02, 98-30



MFR/REFILL #	SIZE	STOCK #	EACH
810007	(5) 60-mL Bottles	40386	\$

For Orion Electrodes (4M KCl sat. with Ag)

Use with: 91-02, 91-03, 91-04, 91-55, 91-56, 91-57, 91-62, 91-63, 91-64, 91-65, 91-66, 91-67, 91-72, 91-80, 96-78

MFR/REFILL #	SIZE	STOCK #	EACH
900011	(5) 60-mL Bottles	40393	\$

For Orion Electrodes (2M KCl sat. with Ag)

Use with: 92-02, 92-03, 92-56, 92-72, 90-04

MFR/REFILL #	SIZE	STOCK #	EACH
900004	(5) 60-mL Bottles	40369	\$



Hi, I'm Abdul K.

“ Our job in Tech Support is to make your job easier. Call or email us for a quick response. ”

We also RECOMMEND...

Electrode Storage Bottles & Solution

Bottle protects standard 12-mm diameter electrodes against bumps and keeps them from drying out between uses. Fill bottle with 50 mL of electrode storage solution.

DESCRIPTION	STOCK #	EACH	QTY 10+
Storage Bottle	40451	\$	\$
Storage Solution, 500 mL	30467		—

